



GTF2B Monoclonal Antibody

Product Code	CSB-MA209230
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant full length protein (316aa) of Human GTF2B.
Raised In	Mouse
Species Reactivity	Human,Mouse,Pig,Rat,Dog
Tested Applications	ELISA,WB,IHC,ICC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	transcription initiation factor IIB S300-II general transcription factor TFIIB RNA polymerase II transcription factor IIB
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	20C10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.