



UHRF1 Monoclonal Antibody

Product Code	CSB-MA439533
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 1aa-408aa of Human UHRF1(N-terminus).
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA, WB, ICC; Recommended dilution: WB:1:500-1:5000, ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%, glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	E3 ubiquitin-protein ligase UHRF1, huNp95, nuclear protein 95, RING finger protein 106, transcription factor ICBP90, nuclear zinc finger protein Np95, inverted CCAAT box-binding protein of 90 kDa
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	10B10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.