



ATIC Monoclonal Antibody

Product Code	CSB-MA136610
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human ATIC protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human,Rabbit
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Relevance	Bifunctional enzyme that catalyzes 2 steps in purine biosynthesis.
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Affinity purified
Isotype	IgG1
Clonality	Monoclonal
Alias	PURH, AICAR, AICARFT, IMPCHASE
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	22G10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.