



SMARCA4 Monoclonal Antibody

Product Code	CSB-MA951350
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human SMARCA4 protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Relevance	Transcriptional coactivator cooperating with nuclear hormone receptors to potentiate transcriptional activation. Component of the CREST-BRG1 complex, a multiprotein complex that regulates promoter activation.
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Affinity purified
Isotype	IgG1
Clonality	Monoclonal
Alias	BRG1, SNF2, SWI2, RTPS2, BAF190, SNF2L4, SNF2LB, hSNF2b, SNF2-BETA
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	4E11-E11-G11