



CRK Monoclonal Antibody

Product Code	CSB-MA289570
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human CRK protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human,Dog,Rabbit,Mouse
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Relevance	The Crk-I and Crk-II forms differ in their biological activities. Crk-II has less transforming activity than Crk-I. Crk-II mediates attachment-induced MAPK8 activation, membrane ruffling and cell motility in a Rac-dependent manner.
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Affinity purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	p38, CRKII
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)