



HSPA2 Monoclonal Antibody

Product Code	CSB-MA054963
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human HSPA1A protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human,Pig,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Relevance	In cooperation with other chaperones, Hsp70s stabilize preexistent proteins against aggregation and mediate the folding of newly translated polypeptides in the cytosol as well as within organelles.
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Affinity purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	HSP72, HSPA1, HSP70I, HSPA1B, HSP70-1, HSP70-1A
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	9D2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.