



ANXA6 Monoclonal Antibody

Product Code	CSB-MA697450
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human ANXA6 protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of ANXA6 and does not cross-react with related proteins.
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Relevance	May associate with CD21. May regulate the release of Ca ²⁺ from intracellular stores.
Form	Purified mouse monoclonal in PBS(pH 7.4) containing with 0.2% sodium azide, 50% glycerol, 0.1% BSA
Purification Method	Affinity purified
Isotype	IgG1
Clonality	Monoclonal
Alias	67 kDa calelectrin; Annexin VI; Annexin-6; Calphobindin-II; Chromobindin-20; Lipocortin VI; Protein III; p68; p70; ANXA6;
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)