



# ANXA6 Monoclonal Antibody

<b>Product Code</b>	CSB-MA697450
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Purified recombinant Human ANXA6 protein fragments expressed in E.coli.
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Specificity</b>	This antibody detects endogenous levels of ANXA6 and does not cross-react with related proteins.
<b>Tested Applications</b>	ELISA,WB;Recommended dilution:WB:1:500-1:5000
<b>Relevance</b>	May associate with CD21. May regulate the release of Ca <sup>2+</sup> from intracellular stores.
<b>Form</b>	Purified mouse monoclonal in PBS(pH 7.4) containing with 0.2% sodium azide, 50% glycerol, 0.1%BSA
<b>Purification Method</b>	Affinity purified
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Alias</b>	67 kDa calelectrin; Annexin VI; Annexin-6; Calphobindin-II; Chromobindin-20; Lipocortin VI; Protein III; p68; p70; ANXA6;
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Clone No.</b>	3D9
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.