



TNFRSF10B Monoclonal Antibody

Product Code	CSB-MA584662
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human TNFRSF10B protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human,Mouse
Specificity	This antibody detects endogenous levels of TNFRSF10B and does not cross-react with related proteins.
Tested Applications	ELISA,WB,ICC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Relevance	Receptor for the cytotoxic ligand TNFSF10/TRAIL.The adapter molecule FADD recruits caspase-8 to the activated receptor.The resulting death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation which initiates the subsequent cascade of caspases (aspartate-specific cysteine proteases) mediating apoptosis. Promotes the activation of NF-kappa-B.Essential for ER stress-induced apoptosis.
Form	Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide,0.1mg/mlBSA and 50% glycerol.
Purification Method	Affinity purified
Isotype	IgG1
Clonality	Monoclonal
Alias	Fas like protein; Apoptosis inducing protein TRICK2A/2B; Apoptosis inducing receptor TRAIL R2; CD 262; CD262; CD262 antigen; Cytotoxic TRAIL receptor 2; Death domain containing receptor for TRAIL/Apo 2L; Death domain containing receptor for TRAIL/Apo2L;
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)