



CDX2 Monoclonal Antibody

Product Code	CSB-MA437110
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human CDX2 protein fragments expressed in E.coli
Raised In	Mouse
Species Reactivity	Human,Monkey,Mouse
Specificity	This antibody detects endogenous levels of CDX2,and does not cross-react with related proteins.
Tested Applications	ELISA,WB,ICC,IP;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200IP:1:200-1:2000
Relevance	Involved in the transcriptional regulation of multiple genes expressed in the intestinal epithelium.Important in broad range of functions from early differentiation to maintenance of the intestinal epithelial lining of both the small and large intestine.
Form	Purified mouse monoclonal in PBS(pH 7.4)containing with 0.2% sodium azide,50% glycerol.
Purification Method	Affinity purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	Caudal type homeo box 2; Caudal type homeo box transcription factor 2; Caudal type homeobox 2; Caudal type homeobox protein 2; caudal type homeobox transcription factor 2; Caudal-type homeobox protein 2; CDX 2; CDX 3; CDX-3; Cdx2; CDX2_HUMAN; CDX3;
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	4B9-B12-E9