



CHD4 Monoclonal Antibody

Product Code	CSB-MA109780
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human CHD4 protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of CHD4 and does not cross-react with related proteins.
Tested Applications	ELISA,WB,IP;Recommended dilution:WB:1:500-1:5000,IP:1:200-1:2000
Relevance	Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones.
Form	Purified mouse monoclonal in PBS(pH 7.4)containing with 0.2% sodium azide,50% glycerol.
Purification Method	Affinity purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	ATP-dependent helicase CHD4; Mi-2 autoantigen 218 kDa protein; Mi2-beta;
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	12C5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.