





RAF1 (Ab-259) Antibody

Product Code	CSB-PA242762
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04049
Immunogen	Peptide sequence around aa.257~261 (S-T-S-T-P) derived from Human RAF.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total Raf1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	C-RAF; C-Raf; CRAF
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAF1

Image



Western blot analysis of extracts from HeLa, K562 and HT29 cells using Raf1 (Ab-259) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Raf1 (Ab-259) Antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.