




RAD51D Antibody


Product Code	CSB-PA019268ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75771
Immunogen	Recombinant Human DNA repair protein RAD51 homolog 4 protein (1-328AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA repair protein RAD51 homolog 4 (R51H3) (RAD51 homolog D) (RAD51-like protein 3) (TRAD), RAD51D, RAD51L3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RAD51D

Image



Western blot

All lanes: DNA repair protein RAD51 homolog 4 antibody at 8µg/ml + Jurkat whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 36, 6, 24, 31, 13, 23, 38 kDa
 Observed band size: 36 kDa



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA019268ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.