




CAPN3 Antibody

Product Code	CSB-PA004497DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20807
Immunogen	Recombinant Human Calpain-3 protein (582-821AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calpain-3 (EC 3.4.22.54) (Calcium-activated neutral proteinase 3) (CANP 3) (Calpain L3) (Calpain p94) (Muscle-specific calcium-activated neutral protease 3) (New calpain 1) (nCL-1), CAPN3, CANP3 CANPL3 NCL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CAPN3

Image



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA004497DSR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.