





IL9 Antibody

Product Code	CSB-PA011675ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15248
Immunogen	Recombinant Human Interleukin-9 protein (19-144AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-9 (IL-9) (Cytokine P40) (T-cell growth factor P40), IL9
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL9

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA011675ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA011675ESR2HU at dilution of 1:100