





C3AR1 Antibody

Product Code	CSB-PA620968ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16581
Immunogen	Recombinant Human C3a anaphylatoxin chemotactic receptor protein (160-340AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	C3a anaphylatoxin chemotactic receptor (C3AR) (C3a-R), C3AR1, AZ3B C3R1 HNFAG09
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	C3AR1

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA620968ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA620968ESR1HU at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.