



POLG2 Antibody

Product Code	CSB-PA883399LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHN1
Immunogen	Recombinant Human DNA polymerase subunit gamma-2, mitochondrial protein (1-485AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase subunit gamma-2, mitochondrial (EC 2.7.7.7) (DNA polymerase gamma accessory 55 kDa subunit) (p55) (Mitochondrial DNA polymerase accessory subunit) (MtPolB) (PolG-beta), POLG2, MTPOLB
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.