



RAB27B Antibody

Product Code	CSB-PA019178HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00194
Immunogen	Recombinant Human Ras-related protein Rab-27B protein (2-218AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-27B (C25KG), RAB27B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAB27B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.