





RETN Antibody

Product Code	CSB-PA07079A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HD89
Immunogen	Recombinant Human Resistin protein (19-108AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Resistin (Adipose tissue-specific secretory factor) (ADSF) (C/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein) (Cysteine-rich secreted protein A12-alpha-like 2) (Cysteine-rich secreted protein FIZZ3), RETN, FIZZ3 HXCP1 RSTN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RETN

Image



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA07079A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA07079A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.