



tuf1 Antibody

Product Code	CSB-PA025342LA01AUD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A3M1F6
Immunogen	Recombinant Acinetobacter baumannii Elongation factor Tu protein (1-396AA)
Raised In	Rabbit
Species Reactivity	Acinetobacter baumannii
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongation factor Tu (EF-Tu), tuf1; tuf2
Immunogen Species	Acinetobacter baumannii (strain ATCC 17978 / CIP 53.77 / LMG 1025 / NCDC KC755 / 5377)
Research Area	Others
Target Names	tuf1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.