





CAPZA2 Antibody

Product Code	CSB-PA004512LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47755
Immunogen	Recombinant Human F-actin-capping protein subunit alpha-2 protein (2-286AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-actin-capping protein subunit alpha-2 (CapZ alpha-2), CAPZA2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CAPZA2

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA004512LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA004512LA01HU at dilution of 1:100