



UGDH Antibody

Product Code	CSB-PA025567HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60701
Immunogen	Recombinant Human UDP-glucose 6-dehydrogenase protein (1-494AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-glucose 6-dehydrogenase (UDP-Glc dehydrogenase) (UDP-GlcDH) (UDPGDH) (EC 1.1.1.22), UGDH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UGDH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.