





AK1 Antibody

Product Code	CSB-PA001508HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00568
Immunogen	Recombinant Human Adenylate kinase isoenzyme 1 protein (1-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenylate kinase isoenzyme 1 (AK 1) (EC 2.7.4.3) (EC 2.7.4.6) (ATP-AMP transphosphorylase 1) (ATP:AMP phosphotransferase) (Adenylate monophosphate kinase) (Myokinase), AK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AK1

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA001508HA01HU at dilution of 1:100



Western Blot

Positive WB detected in: U87 whole cell lysate

All lanes: AK1 antibody at 2.7μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 22 kDa

Observed band size: 22 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.