



cmk Antibody

Product Code	CSB-PA363865HA01EVN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A6I0
Immunogen	Recombinant Escherichia coli Cytidylate kinase protein (1-227AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytidylate kinase (CK) (EC 2.7.4.25) (Cytidine monophosphate kinase) (CMP kinase) (Protein MssA) (p25), cmk, mssA ycaF ycaG
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	cmk
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.