





CCDC36 Antibody

Product Code	CSB-PA816910LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYA8
Immunogen	Recombinant Human Interactor of HORMAD1 protein 1 protein (11-592AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interactor of HORMAD1 protein 1 (Cancer/testis antigen 74) (CT74) (Coiled-coil domain-containing protein 36), CCDC36, IHO1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CCDC36

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA816910LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA816910LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.