



RS1 Antibody

Product Code	CSB-PA020534EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15537
Immunogen	Recombinant Human Retinoschisin protein (24-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinoschisin (X-linked juvenile retinoschisis protein), RS1, XLRS1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.