



CINP Antibody

Product Code	CSB-PA861151LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BW66
Immunogen	Recombinant Human Cyclin-dependent kinase 2-interacting protein (1-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-dependent kinase 2-interacting protein (CDK2-interacting protein), CINP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CINP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.