





NIP7 Antibody

Product Code	CSB-PA897077HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y221
Immunogen	Recombinant Human 60S ribosome subunit biogenesis protein NIP7 homolog protein (1-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosome subunit biogenesis protein NIP7 homolog (KD93), NIP7
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NIP7

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA897077HA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA897077HA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

