



ARSA Antibody

Product Code	CSB-PA002141ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15289
Immunogen	Recombinant Human Arylsulfatase A protein (208-507AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arylsulfatase A (ASA) (EC 3.1.6.8) (Cerebroside-sulfatase) [Cleaved into: Arylsulfatase A component B; Arylsulfatase A component C], ARSA
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ARSA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.