



Agxt Antibody

Product Code	CSB-PA001468LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35423
Immunogen	Recombinant Mouse Serine--pyruvate aminotransferase, mitochondrial protein (25-414AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine--pyruvate aminotransferase, mitochondrial (SPT) (EC 2.6.1.51) (Alanine-glyoxylate aminotransferase) (AGT) (EC 2.6.1.44), Agxt, Agxt1
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Agxt