





BAMBI Antibody

Product Code	CSB-PA615666ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13145
Immunogen	Recombinant Human BMP and activin membrane-bound inhibitor homolog protein (21-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	BMP and activin membrane-bound inhibitor homolog (Non-metastatic gene A protein) (Putative transmembrane protein NMA), BAMBI, NMA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BAMBI


Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA615666ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA615666ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA615666ESR2HU at dilution of 1:100