



GAP43 Antibody

Product Code	CSB-PA009231DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17677
Immunogen	Recombinant Human Neuromodulin protein (1-238AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuromodulin (Axonal membrane protein GAP-43) (Growth-associated protein 43) (Neural phosphoprotein B-50) (pp46), GAP43
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GAP43
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.