



CTNNBIP1 Antibody

Product Code	CSB-PA873645ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NSA3
Immunogen	Recombinant Human Beta-catenin-interacting protein 1 protein (1-81AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-catenin-interacting protein 1 (Inhibitor of beta-catenin and Tcf-4), CTNNBIP1, ICAT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CTNNBIP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.