





SERPINA10 Antibody

Product Code	CSB-PA891944ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK55
Immunogen	Recombinant Human Protein Z-dependent protease inhibitor protein (165-444AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Z-dependent protease inhibitor (PZ-dependent protease inhibitor) (PZI) (SerpA10), SERPINA10, ZPI
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPINA10

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA891944ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA891944ESR1HU at dilution of 1:100