




STAT4 Antibody

Product Code	CSB-PA022813ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14765
Immunogen	Recombinant Human Signal transducer and activator of transcription 4 protein (569-748AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal transducer and activator of transcription 4, STAT4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STAT4

Image



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA022813ESR1HU at dilution of 1:100