




STAT4 Antibody

Product Code	CSB-PA022813ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14765
Immunogen	Recombinant Human Signal transducer and activator of transcription 4 protein (569-748AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal transducer and activator of transcription 4, STAT4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STAT4

Image



Western blot


All lanes: Signal transducer and activator of transcription 4 antibody at 7µg/ml + Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 86 kDa

Observed band size: 86 kDa



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA022813ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.