




SERPINB9 Antibody

Product Code	CSB-PA021077ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50453
Immunogen	Recombinant Human Serpin B9 protein (1-376AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serpin B9 (Cytoplasmic antiproteinase 3) (CAP-3) (CAP3) (Peptidase inhibitor 9) (PI-9), SERPINB9, PI9
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SERPINB9

Image



Western blot

All lanes: Serpin B9 antibody at 3µg/ml + K562 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 42 kDa
 Observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.