





UTS2 Antibody

Product Code	CSB-PA025771ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95399
Immunogen	Recombinant Human Urotensin-2 protein (35-124AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Urotensin-2 (Urotensin II) (U-II) (UII), UTS2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UTS2

Image



Western blot
 All lanes: Urotensin-2 antibody at 8µg/ml + Rat brain tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 15, 17 kDa
 Observed band size: 15 kDa



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA025771ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.