




PARP3 Antibody

Product Code	CSB-PA017467ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6F1
Immunogen	Recombinant Human Poly [ADP-ribose] polymerase 3 protein (294-533AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Poly [ADP-ribose] polymerase 3 (PARP-3) (hPARP-3) (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 3) (ARTD3) (IRT1) (NAD(+) ADP-ribosyltransferase 3) (ADPRT-3) (Poly[ADP-ribose] synthase 3) (pADPRT-3), PARP3, ADPRT3 ADPRTL3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PARP3


Image




Western blot

All lanes: Poly [ADP-ribose] polymerase 3 antibody at 2µg/ml + Mouse heart tissue
Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 60 kDa
Observed band size: 60 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA017467ESR1HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA017467ESR1HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.