




ANXA6 Antibody

Product Code	CSB-PA001847ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08133
Immunogen	Recombinant Human Annexin A6 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Annexin A6 (67 kDa calelectrin) (Annexin VI) (Annexin-6) (Calphobindin-II) (CPB-II) (Chromobindin-20) (Lipocortin VI) (Protein III) (p68) (p70), ANXA6, ANX6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ANXA6

Image



Western blot

All lanes: Annexin A6 antibody at 4µg/ml

Lane 1: Jurkat whole cell lysate

Lane 2: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 76, 73 kDa

Observed band size: 76 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.