




CTBP1 Antibody

Product Code	CSB-PA618772ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13363
Immunogen	Recombinant Human C-terminal-binding protein 1 protein (20-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-terminal-binding protein 1 (CtBP1) (EC 1.1.1.-), CTBP1, CTBP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CTBP1

Image



Western blot
 All lanes: CTBP1 antibody at 6µg/ml + A549 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 48, 47 kDa
 Observed band size: 48 kDa