




ABHD5 Antibody

Product Code	CSB-PA819884ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WTS1
Immunogen	Recombinant Human 1-acylglycerol-3-phosphate O-acyltransferase ABHD5 protein (1-349AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	1-acylglycerol-3-phosphate O-acyltransferase ABHD5 (EC 2.3.1.51) (Abhydrolase domain-containing protein 5) (Lipid droplet-binding protein CGI-58), ABHD5, NCIE2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ABHD5

Image



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA819884ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.