




ACAT2 Antibody

Product Code	CSB-PA861152ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWD1
Immunogen	Recombinant Human Acetyl-CoA acetyltransferase, cytosolic protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acetyl-CoA acetyltransferase, cytosolic (EC 2.3.1.9) (Acetyl-CoA transferase-like protein) (Cytosolic acetoacetyl-CoA thiolase), ACAT2, ACTL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACAT2

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA861152ESR1HU at dilution of 1:100