





AKAP5 Antibody

Product Code	CSB-PA001527ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24588
Immunogen	Recombinant Human A-kinase anchor protein 5 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	A-kinase anchor protein 5 (AKAP-5) (A-kinase anchor protein 79 kDa) (AKAP 79) (H21) (cAMP-dependent protein kinase regulatory subunit II high affinity-binding protein), AKAP5, AKAP79
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	AKAP5

Image



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA001527ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA001527ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.