



ARHGEF9 Antibody

Product Code	CSB-PA002057ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43307
Immunogen	Recombinant Human Rho guanine nucleotide exchange factor 9 protein (307-516AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho guanine nucleotide exchange factor 9 (Collybistin) (PEM-2 homolog) (Rac/Cdc42 guanine nucleotide exchange factor 9), ARHGEF9, ARHDH9 KIAA0424
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGEF9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.