



ARHGEF9 Antibody

Product Code	CSB-PA002057ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43307
Immunogen	Recombinant Human Rho guanine nucleotide exchange factor 9 protein (307-516AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho guanine nucleotide exchange factor 9 (Collybistin) (PEM-2 homolog) (Rac/Cdc42 guanine nucleotide exchange factor 9), ARHGEF9, ARHDH9 KIAA0424
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGEF9