




ARNT2 Antibody

Product Code	CSB-PA864024ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBZ2
Immunogen	Recombinant Human Aryl hydrocarbon receptor nuclear translocator 2 protein (440-690AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aryl hydrocarbon receptor nuclear translocator 2 (ARNT protein 2) (Class E basic helix-loop-helix protein 1) (bHLHe1), ARNT2, BHLHE1 KIAA0307
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ARNT2

Image



Western blot


All lanes: Aryl hydrocarbon receptor nuclear translocator 2 antibody at 4µg/ml + HeLa whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 79, 78 kDa

Observed band size: 79 kDa



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA864024ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.