




AKR1C4 Antibody

Product Code	CSB-PA001545ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17516
Immunogen	Recombinant Human Aldo-keto reductase family 1 member C4 protein (1-323AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aldo-keto reductase family 1 member C4 (EC 1.1.1.-) (3-alpha-HSD1) (3-alpha-hydroxysteroid dehydrogenase type I) (EC 1.1.1.357) (Chlordecone reductase) (CDR) (EC 1.1.1.225) (Dihydrodiol dehydrogenase 4) (DD-4) (DD4) (HAKRA), AKR1C4, CHDR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AKR1C4

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA001545ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.