




# AKR1C4 Antibody

<b>Product Code</b>	CSB-PA001545ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P17516
<b>Immunogen</b>	Recombinant Human Aldo-keto reductase family 1 member C4 protein (1-323AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Aldo-keto reductase family 1 member C4 (EC 1.1.1.-) (3-alpha-HSD1) (3-alpha-hydroxysteroid dehydrogenase type I) (EC 1.1.1.357) (Chlordecone reductase) (CDR) (EC 1.1.1.225) (Dihydrodiol dehydrogenase 4) (DD-4) (DD4) (HAKRA), AKR1C4, CHDR
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	AKR1C4

## Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA001545ESR1HU at dilution of 1:100