



ACP1 Antibody

Product Code	CSB-PA001176ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24666
Immunogen	Recombinant Human Low molecular weight phosphotyrosine protein phosphatase protein (1-158AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low molecular weight phosphotyrosine protein phosphatase (LMW-PTP) (LMW-PTPase) (EC 3.1.3.48) (Adipocyte acid phosphatase) (Low molecular weight cytosolic acid phosphatase) (EC 3.1.3.2) (Red cell acid phosphatase 1), ACP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ACP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.